
ALICA AC Pack Documentation

Release 0.2.0

Marcel Stefko, Kyle M. Douglass

Jun 06, 2018

Contents

1 Javadoc	1
2 Relationship with ALICA and SASS	37
3 Authors	39
4 See Also	41
5 Indices and tables	43

CHAPTER 1

Javadoc

1.1 ch.epfl.leb.alica.acpack.analyzers.autolase

1.1.1 AutoLase

public class **AutoLase** implements Analyzer
Wrapper for Thomas Pengo's implementation of AutoLase algorithm.

Author Marcel Stefko

Constructors

AutoLase

public **AutoLase** (int *threshold*)
Initializes AutoLase with default threshold (120) and averaging (30) values.

Parameters

- **threshold** –

Methods

dispose

public void **dispose** ()

getBatchOutput

public double **getBatchOutput** ()

getIntermittentOutput

```
public double getIntermittentOutput ()
```

getName

```
public String getName ()
```

getShortReturnDescription

```
public String getShortReturnDescription ()
```

getStatusPanel

```
public AnalyzerStatusPanel getStatusPanel ()
```

processImage

```
public void processImage (Object image, int image_width, int image_height, double pixel_size_um, long time_ms)
```

setROI

```
public void setROI (Roi roi)
```

1.1.2 AutoLaseAnalyzer

class **AutoLaseAnalyzer**

This class estimates the density of activations. The density at a particular point relates to the maximum time a certain pixel is “on”, or above a certain threshold. The density is calculated as a moving average 30 frames. The code only works for 2 bytes per pixel cameras for now.

Author Thomas Pengo

Fields

currentDensity

```
double currentDensity
```

running

```
boolean running
```

stopping

boolean **stopping**

Constructors

AutoLaseAnalyzer

public **AutoLaseAnalyzer** (int *threshold*)

Methods

getCurrentValue

public double **getCurrentValue** ()

Returns error signal value from AutoLase

Returns estimated averaged max emitter density

getRawCurrentValue

public double **getRawCurrentValue** ()

Returns raw error signal value from AutoLase

Returns estimated max emitter density for most recent frame

nextImage

public void **nextImage** (ShortProcessor *sp*)

Analyzes next image and adjusts internal state.

Parameters

- **image** – image to be analyzed

setParameters

public void **setParameters** (int *threshold*)

setROI

public void **setROI** (Roi *roi*)

1.1.3 AutoLaseSetupPanel

public class **AutoLaseSetupPanel** extends AnalyzerSetupPanel

Setup panel for AutoLase, allows setup of thresholding and averaging.

Author Marcel Stefko

Constructors

AutoLaseSetupPanel

```
public AutoLaseSetupPanel()  
    Creates new form SetupPanel
```

Methods

getName

```
public String getName()
```

initAnalyzer

```
public Analyzer initAnalyzer()
```

toString

```
public String toString()
```

1.2 ch.epfl.leb.alica.acpack.analyzers.defcon

1.2.1 DEFCoNSetupPanel

```
public class DEFCoNSetupPanel extends AnalyzerSetupPanel  
    The DEFCoN analyzer setup panel.
```

Author Kyle M. Douglass

Constructors

DEFCoNSetupPanel

```
public DEFCoNSetupPanel()  
    Creates new form SetupPanel
```

Methods

getName

```
public String getName()
```

initAnalyzer

```
public Analyzer initAnalyzer()  
    Initializes the analyzer using the properties from the setup panel.
```

1.2.2 Defcon

```
public class Defcon implements Analyzer  
    A fluorescent spot counter derived from the DEFCoN package.
```

Author Kyle M. Douglass

See also: [DEFCoN-ImageJ](#)

Constructors

Defcon

```
public Defcon (String pathToModel)  
    Initializes the DEFCoN analyzer.
```

Parameters

- **pathToModel** – The path to the DEFCoN network model.

Methods

dispose

```
public void dispose()  
    Cleans up the analyzer when it's finished.
```

finalize

```
protected void finalize()  
    Failsafe in case the predictor has not been closed at the point of garbage collection.
```

Throws

- [java.lang.Throwable](#) –

getBatchOutput

```
public double getBatchOutput()  
    Returns the averaged DEFCoN count value since the last call. Double.NaN is returned if there is no new count since the previous call.
```

Returns The averaged DEFCoN count.

getBoxSize

```
public int getBoxSize()  
    Returns the current square kernel size for maximum local counts.
```

Returns The kernel size for computing the maximum local count.

getIntermittentOutput

```
public double getIntermittentOutput()  
    Returns the intermittent output of the analyzer.
```

Returns The analyzer's current output value.

getName

```
public String getName()  
    Returns the name of the DEFCoN analyzer.
```

Returns The name of the DEFCoN analyzer.

getShortReturnDescription

```
public String getShortReturnDescription()  
    Returns a short description of the values returned by the DEFCoN analyzer.
```

Returns A short description of the values returned by the DEFCoN analyzer.

getStatusPanel

```
public AnalyzerStatusPanel getStatusPanel()  
    Returns the analyzer's status panel that will be displayed in the GUI. If no panel is implemented, this method  
    should return null. In this case, the corresponding space in the MonitorGUI will appear blank.
```

Returns The status panel of the DEFCoN analyzer or null.

isLiveModeOn

```
public boolean isLiveModeOn()  
    True if live mode is on, false otherwise.
```

Returns True if live mode is on, false otherwise

isMaxLocalCount

```
public boolean isMaxLocalCount()  
    True if the analyzer is computing the maximum local count.
```

Returns True if the analyzer is computing the maximum local count.

liveModeOff

```
public void liveModeOff()  
    Turns off the live view of the density map.
```

liveModeOn

```
public void liveModeOn()  
    Turns on the live view of the density map.
```

maxLocalCountOff

```
public void maxLocalCountOff()  
    Turns off the live view of the density map.
```

maxLocalCountOn

```
public void maxLocalCountOn()  
    Turns on the live view of the density map.
```

processImage

```
public void processImage (Object image, int width, int height, double pixelSizeUm, long timeMs)  
Processes an image and adjusts the analyzer's internal state to reflect the results of the calculation. This method  
is called after each new image acquisition by the AnalysisWorker. You can use the synchronized(this) statement  
within the body of an implementation of an Analyzer to ensure that no output readout happens during code  
execution.
```

Parameters

- **image** – The image to be processed as 1D raw pixel data.
- **width** – Image width in pixels.
- **height** – Image height in pixels.
- **pixelSizeUm** – Length of a side of a square pixel in micrometers.
- **timeMs** – Image acquisition time in milliseconds.

setBoxSize

```
public void setBoxSize (int boxSize)  
Sets the square kernel size for computing the maximum local count.
```

Parameters

- **boxSize** – The kernel size for computing the maximum local count.

setROI

```
public void setROI (Roi roi)
```

updateLiveView

```
public void updateLiveView()  
    Updates the live viewer.
```

1.2.3 DefconStatusPanel

```
public class DefconStatusPanel extends AnalyzerStatusPanel  
    Status panel of the DEFCoN analyzer.
```

Author Kyle M. Douglass

Constructors

DefconStatusPanel

```
public DefconStatusPanel (Defcon defcon)  
    Creates new form DefconStatusPanel.
```

1.2.4 DefconTest

```
public class DefconTest  
    Unit tests for the Defcon class.
```

Author Kyle M. Douglass

Methods

setUp

```
public void setUp()
```

testGetSetBoxSize

```
public void testGetSetBoxSize()  
    Gets/sets the box size field for the maximum local count.
```

testMaxLocalCount

```
public void testMaxLocalCount()  
    Toggles the maximum local count feature.
```

testNegativeBoxSize

```
public void testNegativeBoxSize()  
    Ensures that the boxSize value is greater than 1.
```

testOddBoxSize

```
public void testOddBoxSize()  
    Ensures that the boxSize value is odd.
```

testProcessImage

```
public void testProcessImage()  
    Ensures that the processImage() method is called without errors.
```

testProcessMaxLocalCount

```
public void testProcessMaxLocalCount()  
    Analyzer predicts the maximum local count instead of the density map.
```

1.3 ch.epfl.leb.alica.acpack.analyzers.integrator

1.3.1 Integrator

```
public class Integrator implements Analyzer  
    Analyzer which outputs the average pixel value per frame. The average is taken over the area of the image (or  
    ROI) in units of squared pixels.
```

Author Marcel Stefko

Constructors**Integrator**

```
public Integrator()
```

Methods**dispose**

```
public void dispose()
```

getBatchOutput

```
public double getBatchOutput()
```

getIntermittentOutput

```
public double getIntermittentOutput()
```

getName

```
public String getName ()
```

getShortReturnDescription

```
public String getShortReturnDescription ()
```

getStatusPanel

```
public AnalyzerStatusPanel getStatusPanel ()
```

processImage

```
public void processImage (Object image, int image_width, int image_height, double pixel_size_um, long time_ms)
```

Computes the average of the pixel values taken over the image (or ROI).

Parameters

- `image` –
- `image_width` –
- `image_height` –
- `pixel_size_um` –
- `time_ms` –

setROI

```
public void setROI (Roi roi)
```

1.3.2 IntegratorSetupPanel

```
public class IntegratorSetupPanel extends AnalyzerSetupPanel
```

Empty panel

Author stefko

Constructors

IntegratorSetupPanel

```
public IntegratorSetupPanel ()
```

Creates new form IntegratorSetupPanel

Methods

getName

```
public String getName()
```

initAnalyzer

```
public Analyzer initAnalyzer()
```

1.3.3 IntegratorTest

```
public class IntegratorTest
```

Author douglass

Methods

setUp

```
public void setUp()
```

testGetBatchOutput

```
public void testGetBatchOutput()
```

Test of getBatchOutput method, of class Integrator.

testProcessImage

```
public void testProcessImage()
```

Test of processImage method, of class Integrator.

1.4 ch.epfl.leb.alica.acpack.analyzers.quickpalm

1.4.1 MyDialogs

```
class MyDialogs
```

Fields

attach

```
boolean attach
```

buffer

int **buffer**

cal_z

double **cal_z**

calfile

java.lang.String **calfile**

fwHM

double **fwHM**

height

int **height**

imagedir

java.lang.String **imagedir**

imp

ImagePlus **imp**

imtitle

String **imtitle**

is3d

boolean **is3d**

magn

double **magn**

maxpart

int **maxpart**

maxsizedouble **maxsize****minsize**double **minsize****model**java.lang.String **model****models**java.lang.String[] **models****nimchars**int **nimchars****nimstart**int **nimstart****nrois**int **nrois****nslices**int **nslices****part_divergence**boolean **part_divergence****part_extrainfo**boolean **part_extrainfo****pattern**java.lang.String **pattern**

pixelsize

```
double pixelsize
```

prefix

```
java.lang.String prefix
```

prefs

```
ij.Prefs prefs
```

ptablefile

```
java.lang.String ptablefile
```

pthrsh

```
double pthrsh
```

rmanager

```
RoiManager rmanager
```

rois

```
Roi[] rois
```

saturation

```
double saturation
```

smartsnr

```
boolean smartsnr
```

snr

```
int snr
```

sufix

```
java.lang.String sufix
```

symmetry

double **symmetry**

threads

int **threads**

view

boolean **view**

view_mode

java.lang.String **view_mode**

view_modes

java.lang.String[] **view_modes**

viewer_accumulate

int **viewer_accumulate**

viewer_do3d

boolean **viewer_do3d**

viewer_doConvolve

boolean **viewer_doConvolve**

viewer_doMovie

boolean **viewer_doMovie**

viewer_fwhm

double **viewer_fwhm**

viewer_is8bit

boolean **viewer_is8bit**

viewer_mergeabove

```
double viewer_mergeabove
```

viewer_mergebelow

```
double viewer_mergebelow
```

viewer_oheight

```
int viewer_oheight
```

viewer_owidth

```
int viewer_owidth
```

viewer_tpixelsize

```
double viewer_tpixelsize
```

viewer_update

```
int viewer_update
```

viewer_zstep

```
double viewer_zstep
```

waittime

```
int waittime
```

width

```
int width
```

window

```
int window
```

Methods

analyseParticles

```
public boolean analyseParticles (MyFunctions f)
```

1.4.2 MyFunctions

```
class MyFunctions
```

Fields

cal3d_center

```
int cal3d_center
```

cal3d_wmh

```
double[] cal3d_wmh
```

cal3d_z

```
double[] cal3d_z
```

caltable

```
ResultsTable caltable
```

debug

```
boolean debug
```

dtable

```
ResultsTable dtable
```

gblur

```
GaussianBlur gblur
```

live_view

```
ImagePlus live_view
```

ptable

ResultsTable **ptable**

ptable_lock

java.util.concurrent.locks.Lock **ptable_lock**

Constructors

MyFunctions

public **MyFunctions** (boolean *live_view*)

Methods

argmax

int **argmax** (double[] *arr*)

argmin

int **argmin** (double[] *arr*)

clearRegion

void **clearRegion** (double *thrsh*, ImageProcessor *ip*, boolean[][] *mask*, int *xstart*, int *xend*, int *ystart*, int *yend*)

detectParticles

int **detectParticles** (ImageProcessor *ip*, *MyDialogs* *dg*, int *nframe*)

Particle finding method, will search the image for particles.

Parameters

- **ip** – image to search for particles on
- **dg** – dialog manager
- **nframe** – the frame index corresponding to this image

dispose

public void **dispose** ()

getClosest

```
int getClosest (double value, double[] arr, int center)
```

getMaxPositions

```
int[] getMaxPositions (ImageProcessor ip)
```

getNextImage

```
ImagePlus getNextImage (MyDialogs dg, int frame)
```

Grabs a new image on an observed folder on the case of the analysis being attached to the acquisition.

Parameters

- **dg** – dialog manager
- **frame** – frame index to search for on the folder

Returns found image

getParticle

```
boolean getParticle (ImageProcessor ip, boolean[][] mask, int[] maxs, MyDialogs dg, ResultsTable ptable,  
int nframe)
```

Particle analysis method, called for each particle candidate found by detectParticles.

Parameters

- **ip** – image to search for particles on
- **dg** – dialog manager
- **nframe** – the frame index corresponding to this image

getParticleForCalibration

```
double[] getParticleForCalibration (ImageProcessor ip, MyDialogs dg, int xstart, int xend, int ystart, int yend)
```

getZ

```
double getZ (double wmh)
```

Given a calculated width-minus-height (wmh) converts this value into the corresponding coordinate in Z by comparing against the loaded Z-calibration table. Only used if the particle is disturbed by astigmatism.

Parameters

- **wmh** – the width-minus-height of a particle

Returns corresponding z-position value, will return 9999 if wmh is out of limits

initialize3d

```
void initialize3d()  
    Loads values from the calibration table into the cal3d_* arrays.
```

log

```
void log (java.lang.String txt)
```

mean

```
double mean (double[] array, int start, int stop)
```

movingMean

```
double[] movingMean (double[] arr, int window)
```

showTable

```
void showTable ()
```

1.4.3 QuickPalm

public class **QuickPalm** implements Analyzer
Produces a localization count per area using QuickPALM.

Author Marcel Stefko

Constructors

QuickPalm

```
public QuickPalm (boolean live_view)  
Implementation of the QuickPALM algorithm as an analyzer, which produces particle count as output.
```

Parameters

- **live_view** – if true, live view of particle positions is shown

Methods

dispose

```
public void dispose ()
```

getBatchOutput

```
public double getBatchOutput ()
```

getIntermittentOutput

```
public double getIntermittentOutput ()
```

getName

```
public String getName ()
```

getShortReturnDescription

```
public String getShortReturnDescription ()
```

getStatusPanel

```
public AnalyzerStatusPanel getStatusPanel ()
```

processImage

```
public void processImage (Object image, int image_width, int image_height, double pixel_size_um, long time_ms)
```

setROI

```
public void setROI (Roi roi)
```

1.4.4 QuickPalmCore**public class QuickPalmCore**

Wrapper for QuickPalm ImageJ plugin functionality developed by Ricardo Henriques @ Instituto de Medicina Molecular (PT) / Institut Pasteur (FR).

Author Marcel Stefko

Constructors**QuickPalmCore**

```
public QuickPalmCore (boolean live_view)
```

Initializes the core and launches a dialog for parameter setup.

Parameters

- **live_view** – if true, live view of particle positions will be shown

Methods

dispose

```
public void dispose ()  
    Close preview window if opened.
```

processImage

```
public int processImage (ImageProcessor ip, int frame)  
    Counts particles in the image.
```

Parameters

- **ip** – image to be processed
- **frame** – id of the image

Returns no. of detected particles

1.4.5 QuickPalmSetupPanel

```
public class QuickPalmSetupPanel extends AnalyzerSetupPanel
```

Author stefko

Constructors

QuickPalmSetupPanel

```
public QuickPalmSetupPanel ()  
    Creates new form QuickPalmSetupPanel
```

Methods

getName

```
public String getName ()
```

initAnalyzer

```
public Analyzer initAnalyzer ()
```

1.5 ch.epfl.leb.alica.acpack.analyzers.spotcounter

1.5.1 FindLocalMaxima

```
public class FindLocalMaxima
```

Find local maxima in an Image (or ROI) using the algorithm described in Neubeck and Van Gool. Efficient non-

maximum suppression. Pattern Recognition (2006) vol. 3 pp. 850-855 Jonas Ries brought this to my attention and send me C code implementing one of the described algorithms

Methods

FindMax

public static **Polygon FindMax** (ImageProcessor *iProc*, Roi *roi*, int *n*, int *threshold*, *FilterType filterType*)
 Static utility function to find local maxima in an Image

Parameters

- **iProc** –
 - ImageProcessor object in which to look for local maxima
- **roi** –
 - region of interest to which the analysis is constrained
- **n** –
 - minimum distance to other local maximum
- **threshold** –
 - value below which a maximum will be rejected
- **filterType** –
 - Prefilter the image. Either none or Gaussian1_5

Returns Polygon with maxima

noiseFilter

public static **Polygon noiseFilter** (ImageProcessor *iProc*, **Polygon inputPoints**, int *threshold*)

Parameters

- **iProc** –
- **inputPoints** –
- **threshold** –

1.5.2 FindLocalMaxima.FilterType

public enum **FilterType**

Different filters for image preprocessing.

Enum Constants

GAUSSIAN1_5

public static final *FindLocalMaxima.FilterType GAUSSIAN1_5*

NONE

```
public static final FindLocalMaxima.FilterType NONE  
    No preprocessing.
```

1.5.3 SpotCounter

```
public class SpotCounter implements Analyzer
```

Author stefko

Constructors

SpotCounter

```
public SpotCounter (int noise_tolerance, int box_size, boolean live_view)  
    Initialize the analyzer
```

Parameters

- **noise_tolerance** – required height of peak around surroundings
- **box_size** – size of the scanning box in pixels
- **live_view** – if true, live preview is shown

Methods

dispose

```
public void dispose ()
```

getBatchOutput

```
public double getBatchOutput ()
```

getIntermittentOutput

```
public double getIntermittentOutput ()
```

getName

```
public String getName ()
```

getShortReturnDescription

```
public String getShortReturnDescription ()
```

getStatusPanel

```
public AnalyzerStatusPanel getStatusPanel ()
```

processImage

```
public void processImage (Object image, int image_width, int image_height, double pixel_size_um, long time_ms)
```

setROI

```
public void setROI (Roi roi)
```

1.5.4 SpotCounterCore

```
public class SpotCounterCore  
Core of the SpotCounter algorithm
```

Author Nico Stuurman

Constructors**SpotCounterCore**

```
public SpotCounterCore (int noiseTolerance, int boxSize, boolean live_mode)
```

Parameters

- **noiseTolerance** – minimum peak value
- **boxSize** – size of scanning box
- **live_mode** – if true, live preview is shown

Methods**analyze**

```
public HashMap<String, Double> analyze (ImageProcessor ip)  
Analyzes the image and returns information about current state.
```

Parameters

- **ip** –
 - image to be analyzed

Returns ResultsTable which contains information about analysis results.

dispose

```
public void dispose ()  
Hide live view window if it exists.
```

getBoxSize

```
public int getBoxSize ()
```

getNoiseTolerance

```
public int getNoiseTolerance ()
```

isLiveModeOn

```
public boolean isLiveModeOn ()
```

Returns true if live mode is on, false otherwise

liveModeOff

```
public void liveModeOff ()
```

Turns off live viewing of SpotCounter analysis.

liveModeOn

```
public void liveModeOn ()
```

Turns on live viewing of SpotCounter analysis.

setParams

```
public void setParams (int noiseTolerance, int boxSize)
```

Set new parameters for the analysis

Parameters

- **noiseTolerance** – minimum peak value
- **boxSize** – size of scanning box

setROI

```
public void setROI (Roi roi)
```

Constrain analysis to given ROI.

Parameters

- **roi** – ROI to constrain analysis to

1.5.5 SpotCounterSetupPanel

```
public class SpotCounterSetupPanel extends AnalyzerSetupPanel
```

Setup panel to initialize the SpotCounter

Author Marcel Stefko

Constructors

SpotCounterSetupPanel

```
public SpotCounterSetupPanel()  
    Creates new form SetupPanel
```

Methods

getName

```
public String getName()
```

initAnalyzer

```
public Analyzer initAnalyzer()
```

1.5.6 SpotCounterStatusPanel

```
public class SpotCounterStatusPanel extends AnalyzerStatusPanel  
    Status panel of SpotCounter, enables modification of threshold in real time.
```

Author Marcel Stefko

Constructors

SpotCounterStatusPanel

```
public SpotCounterStatusPanel (SpotCounterCore core)  
    Creates new form SpotCounterStatusPanel
```

Parameters

- **core** – SpotCounterCore

1.6 ch.epfl.leb.alica.acpack.controllers.inverter

1.6.1 InvertController

```
public class InvertController implements Controller  
    Controller which inverts and scales the input using 1/x function. (high input -> low output, low input -> high output)
```

Author Marcel Stefko

Fields

maximum

protected double **maximum**
Maximal possible output value.

Constructors

InvertController

public **InvertController** (double *maximum*, double *value_at_1_mw*)
Initializes the InvertController

Parameters

- **maximum** – max output value
- **value_at_1_mw** – what is the value of input that you want to cause an output value of 1.0 (scaling constant)

Methods

getCurrentOutput

public double **getCurrentOutput** ()

getName

public String **getName** ()

getSetpoint

public double **getSetpoint** ()

getStatusPanel

public ControllerStatusPanel **getStatusPanel** ()

nextValue

public double **nextValue** (double *value*)

setSetpoint

```
public void setSetpoint (double value)
```

Sets the scaling constant, since it basically fulfills the role of setpoint for this controller.

Parameters

- **value** – new scaling constant

1.6.2 InverterSetupPanel

```
public class InverterSetupPanel extends ControllerSetupPanel
```

Setup panel for the InvertController

Author Marcel Stefko

Constructors

InverterSetupPanel

```
public InverterSetupPanel ()
```

Creates new form InverterSetupPanel

Methods

getName

```
public String getName ()
```

initController

```
public Controller initController (double max_controller_output, double tick_rate_ms)
```

1.7 ch.epfl.leb.alica.acpack.controllers.manual

1.7.1 ManualController

```
public class ManualController implements Controller
```

Manual controller. Output is equal to setpoint value, any input is ignored.

Author Marcel Stefko

Fields

maximum

```
protected double maximum
```

Maximal possible output value.

setpoint

```
protected double setpoint  
    Output value
```

Constructors

ManualController

```
public ManualController (double maximum, double initial_output)  
    Initialize with maximal output value
```

Parameters

- **maximum** – max output value
- **initial_output** – starting value of output

Methods

getCurrentOutput

```
public double getCurrentOutput ()
```

getName

```
public String getName ()
```

getSetpoint

```
public double getSetpoint ()
```

getStatusPanel

```
public ControllerStatusPanel getStatusPanel ()
```

nextValue

```
public double nextValue (double value)
```

setSetpoint

```
public void setSetpoint (double new_setpoint)
```

1.7.2 ManualSetupPanel

```
public class ManualSetupPanel extends ControllerSetupPanel
```

Author stefko

Constructors

ManualSetupPanel

```
public ManualSetupPanel()
```

Creates new form ManualSetupPanel

Methods

getName

```
public String getName()
```

initController

```
public Controller initController(double max_controller_output, double tick_rate_ms)
```

1.8 ch.epfl.leb.alica.acpack.controllers.pi

1.8.1 PI_SetupPanel

```
public class PI_SetupPanel extends ControllerSetupPanel
```

Author stefko

Constructors

PI_SetupPanel

```
public PI_SetupPanel()
```

Creates new form PI_SetupPanel

Methods

getName

```
public String getName()
```

initController

```
public Controller initController(double max_controller_output, double tick_rate_ms)
```

1.8.2 PI_StatusPanel

public class **PI_StatusPanel** extends ControllerStatusPanel

Author stefko

Constructors

PI_StatusPanel

public **PI_StatusPanel** (*PI_controller controller*)

Creates new form PI_StatusPanel

Methods

setValuesDisplay

public void **setValuesDisplay** (double *P*, double *I*)

1.8.3 PI_controller

public class **PI_controller** implements Controller

Simple implementation of PI controller with output constraining and windup prevention.

Author Marcel Stefko

Fields

I

protected double **I**

Integral component

P

protected double **P**

Proportional component

current_output

protected double **current_output**

Last calculated output of the controller

is_blocked

protected boolean **is_blocked**

Constructors

PI_controller

```
public PI_controller (double P, double I_per_second, double max_output, double sampling_period_s)  
    Initialize the PI controller.
```

Parameters

- **P** – proportional component
- **I_per_second** – integral component (per second)
- **max_output** – maximal output value
- **sampling_period_s** – controller tick rate in seconds

Methods

block

```
public void block ()  
    Temporarily stops the controller from taking in input, and forces output to be 0.
```

getCurrentOutput

```
public double getCurrentOutput ()
```

getName

```
public String getName ()
```

getSetpoint

```
public double getSetpoint ()
```

getStatusPanel

```
public ControllerStatusPanel getStatusPanel ()
```

nextValue

```
public double nextValue (double value)
```

setSetpoint

```
public void setSetpoint (double new_setpoint)
```

unblock

```
public void unblock()  
    Resets integral before unblocking the output
```

1.9 ch.epfl.leb.alica.acpack.controllers.selftuningpi

1.9.1 SelfTuningController

```
public class SelfTuningController extends PI_controller
```

A self-tuning implementation of the PI controller. It waits for 10 cycles, and over next 20 cycles generates a step pulse and measures response. From this response it estimates the P and I components of the PID controller.

Author Marcel Stefko

Constructors

SelfTuningController

```
public SelfTuningController (double max_output, double sampling_period_s)
```

Initialize a new SelfTuningController

Parameters

- **max_output** – max output
- **sampling_period_s** – tick rate of controller in seconds
- **step_height** – how big step pulse should be generated
- **p_factor** – scaling factor for P in tuning
- **i_factor** – scaling factor for I in tuning

Methods

getName

```
public String getName ()
```

getStatusPanel

```
public ControllerStatusPanel getStatusPanel ()
```

nextValue

```
public double nextValue (double value)
```

recalibrate

```
public final void recalibrate (double step_height, double p_factor, double i_factor)
```

1.9.2 SelfTuningSetupPanel

```
public class SelfTuningSetupPanel extends ControllerSetupPanel
```

Author stefko

Constructors

SelfTuningSetupPanel

```
public SelfTuningSetupPanel ()
```

Creates new form SelfTuningSetupPanel

Methods

getName

```
public String getName ()
```

initController

```
public Controller initController (double max_controller_output, double tick_rate_ms)
```

1.9.3 SelfTuningStatusPanel

```
public class SelfTuningStatusPanel extends ControllerStatusPanel
```

Author stefko

Constructors

SelfTuningStatusPanel

```
public SelfTuningStatusPanel (SelfTuningController controller)
```

Creates new form SeltTuningStatusPanel

Methods

setCalibrationStatusDisplay

```
public void setCalibrationStatusDisplay (String text)
```


setValuesDisplay

```
public void setValuesDisplay (double P, double I)
```

This package contains the implementations for the analyzers and controllers in ALICA, the automatic illumination control package for light microscopy. It is a run-time dependency of ALICA. These implementations are maintained separately of ALICA because they are used in SASS as well.

CHAPTER 2

Relationship with ALICA and SASS



- Dependency at compile time, parent class files are also copied into final JAR.
- Dependency at compile time, parent class files NOT included in final JAR.
- Dynamically loaded at runtime, not a compile-time dependency.

CHAPTER 3

Authors

- Marcel Štefko
- Kyle M. Douglass

CHAPTER 4

See Also

- [ALICA - Automated Laser Illumination Control Algorithm](#)
- [SASS - SMLM Acquisition Simulation Software](#)

CHAPTER 5

Indices and tables

- genindex

Index

A

analyseParticles(MyFunctions) (Java method), 17
analyze(ImageProcessor) (Java method), 25
argmax(double[]) (Java method), 18
argmin(double[]) (Java method), 18
attach (Java field), 11
AutoLase (Java class), 1
AutoLase(int) (Java constructor), 1
AutoLaseAnalyzer (Java class), 2
AutoLaseAnalyzer(int) (Java constructor), 3
AutoLaseSetupPanel (Java class), 3
AutoLaseSetupPanel() (Java constructor), 4

B

block() (Java method), 33
buffer (Java field), 12

C

cal3d_center (Java field), 17
cal3d_wmh (Java field), 17
cal3d_z (Java field), 17
cal_z (Java field), 12
calfile (Java field), 12
caltable (Java field), 17
ch.epfl.leb.alica.acpack.analyzers.autolase (package), 1
ch.epfl.leb.alica.acpack.analyzers.defcon (package), 4
ch.epfl.leb.alica.acpack.analyzers.integrator (package), 9
ch.epfl.leb.alica.acpack.analyzers.quickpalm (package), 11
ch.epfl.leb.alica.acpack.analyzers.spotcounter (package), 22
ch.epfl.leb.alica.acpack.controllers.inverter (package), 27
ch.epfl.leb.alica.acpack.controllers.manual (package), 29
ch.epfl.leb.alica.acpack.controllers.pi (package), 31
ch.epfl.leb.alica.acpack.controllers.selftuningpi (package), 34
clearRegion(double, ImageProcessor, boolean[][], int, int, int, int) (Java method), 18
current_output (Java field), 32

currentDensity (Java field), 2

D

debug (Java field), 17
Defcon (Java class), 5
Defcon(String) (Java constructor), 5
DEFCoNSetupPanel (Java class), 4
DEFCoNSetupPanel() (Java constructor), 4
DefconStatusPanel (Java class), 8
DefconStatusPanel(Defcon) (Java constructor), 8
DefconTest (Java class), 8
detectParticles(ImageProcessor, MyDialogs, int) (Java method), 18
dispose() (Java method), 1, 5, 9, 18, 20, 22, 24, 25
dtable (Java field), 17

F

FilterType (Java enum), 23
finalize() (Java method), 5
FindLocalMaxima (Java class), 22
FindMax(ImageProcessor, Roi, int, int, FilterType) (Java method), 23
fwhm (Java field), 12

G

GAUSSIAN1_5 (Java field), 23
gblur (Java field), 17
getBatchOutput() (Java method), 1, 5, 9, 21, 24
getBoxSize() (Java method), 6, 26
getClosest(double, double[], int) (Java method), 19
getCurrentOutput() (Java method), 28, 30, 33
getCurrentValue() (Java method), 3
getIntermittentOutput() (Java method), 2, 6, 9, 21, 24
getMaxPositions(ImageProcessor) (Java method), 19
getName() (Java method), 2, 4, 6, 10, 11, 21, 22, 24, 27–31, 33–35
getNextImage(MyDialogs, int) (Java method), 19
getNoiseTolerance() (Java method), 26
getParticle(ImageProcessor, boolean[][], int[], MyDialogs, ResultsTable, int) (Java method), 19

getParticleForCalibration(ImageProcessor, int, int, int, int) (Java method), 19
getRawCurrentValue() (Java method), 3
getSetpoint() (Java method), 28, 30, 33
getShortReturnDescription() (Java method), 2, 6, 10, 21, 24
getStatusPanel() (Java method), 2, 6, 10, 21, 25, 28, 30, 33, 34
getZ(double) (Java method), 19

H

height (Java field), 12

I

I (Java field), 32
imagedir (Java field), 12
imp (Java field), 12
imtitle (Java field), 12
initAnalyzer() (Java method), 4, 5, 11, 22, 27
initController(double, double) (Java method), 29, 31, 35
initialize3d() (Java method), 20
Integrator (Java class), 9
Integrator() (Java constructor), 9
IntegratorSetupPanel (Java class), 10
IntegratorSetupPanel() (Java constructor), 10
IntegratorTest (Java class), 11
InvertController (Java class), 27
InvertController(double, double) (Java constructor), 28
InverterSetupPanel (Java class), 29
InverterSetupPanel() (Java constructor), 29
is3d (Java field), 12
is_blocked (Java field), 32
isLiveModeOn() (Java method), 6, 26
isMaxLocalCount() (Java method), 6

L

live_view (Java field), 17
liveModeOff() (Java method), 7, 26
liveModeOn() (Java method), 7, 26
log(java.lang.String) (Java method), 20

M

magn (Java field), 12
ManualController (Java class), 29
ManualController(double, double) (Java constructor), 30
ManualSetupPanel (Java class), 31
ManualSetupPanel() (Java constructor), 31
maximum (Java field), 28, 29
maxLocalCountOff() (Java method), 7
maxLocalCountOn() (Java method), 7
maxpart (Java field), 12
maxsize (Java field), 13
mean(double[], int, int) (Java method), 20

minsize (Java field), 13
model (Java field), 13
models (Java field), 13
movingMean(double[], int) (Java method), 20
MyDialogs (Java class), 11
MyFunctions (Java class), 17
MyFunctions(boolean) (Java constructor), 18

N

nextImage(ShortProcessor) (Java method), 3
nextValue(double) (Java method), 28, 30, 33, 34
nimchars (Java field), 13
nimstart (Java field), 13
noiseFilter(ImageProcessor, Polygon, int) (Java method), 23
NONE (Java field), 24
nrois (Java field), 13
nslices (Java field), 13

P

P (Java field), 32
part_divergence (Java field), 13
part_extrainfo (Java field), 13
pattern (Java field), 13
PI_controller (Java class), 32
PI_controller(double, double, double, double) (Java constructor), 33
PI_SetupPanel (Java class), 31
PI_SetupPanel() (Java constructor), 31
PI_StatusPanel (Java class), 32
PI_StatusPanel(PI_controller) (Java constructor), 32
pixelsize (Java field), 14
prefix (Java field), 14
prefs (Java field), 14
processImage(ImageProcessor, int) (Java method), 22
processImage(Object, int, int, double, long) (Java method), 2, 7, 10, 21, 25
ptable (Java field), 18
ptable_lock (Java field), 18
ptablefile (Java field), 14
pthrsh (Java field), 14

Q

QuickPalm (Java class), 20
QuickPalm(boolean) (Java constructor), 20
QuickPalmCore (Java class), 21
QuickPalmCore(boolean) (Java constructor), 21
QuickPalmSetupPanel (Java class), 22
QuickPalmSetupPanel() (Java constructor), 22

R

recalibrate(double, double, double) (Java method), 35
rmanager (Java field), 14

rois (Java field), 14
running (Java field), 2

S

saturation (Java field), 14
SelfTuningController (Java class), 34
SelfTuningController(double, double) (Java constructor), 34
SelfTuningSetupPanel (Java class), 35
SelfTuningSetupPanel() (Java constructor), 35
SelfTuningStatusPanel (Java class), 35
SelfTuningStatusPanel(SelfTuningController) (Java constructor), 35
setBoxSize(int) (Java method), 7
setCalibrationStatusDisplay(String) (Java method), 35
setParameters(int) (Java method), 3
setParams(int, int) (Java method), 26
setpoint (Java field), 30
setROI(Roi) (Java method), 2, 3, 7, 10, 21, 25, 26
setSetpoint(double) (Java method), 29, 30, 33
setUp() (Java method), 8, 11
setValuesDisplay(double, double) (Java method), 32, 36
showTable() (Java method), 20
smartsnr (Java field), 14
snr (Java field), 14
SpotCounter (Java class), 24
SpotCounter(int, int, boolean) (Java constructor), 24
SpotCounterCore (Java class), 25
SpotCounterCore(int, int, boolean) (Java constructor), 25
SpotCounterSetupPanel (Java class), 26
SpotCounterSetupPanel() (Java constructor), 27
SpotCounterStatusPanel (Java class), 27
SpotCounterStatusPanel(SpotCounterCore) (Java constructor), 27
stopping (Java field), 3
suffix (Java field), 14
symmetry (Java field), 15

T

testGetBatchOutput() (Java method), 11
testGetBoxSize() (Java method), 8
testMaxLocalCount() (Java method), 8
testNegativeBoxSize() (Java method), 8
testOddBoxSize() (Java method), 9
testProcessImage() (Java method), 9, 11
testProcessMaxLocalCount() (Java method), 9
threads (Java field), 15
toString() (Java method), 4

U

unblock() (Java method), 34
updateLiveView() (Java method), 8

V

view (Java field), 15
view_mode (Java field), 15
view_modes (Java field), 15
viewer_accumulate (Java field), 15
viewer_do3d (Java field), 15
viewer_doConvolve (Java field), 15
viewer_doMovie (Java field), 15
viewer_fwhm (Java field), 15
viewer_is8bit (Java field), 15
viewer_mergeabove (Java field), 16
viewer_mergebelow (Java field), 16
viewer_oheight (Java field), 16
viewer_owidth (Java field), 16
viewer_tpixelsize (Java field), 16
viewer_update (Java field), 16
viewer_zstep (Java field), 16

W

waittime (Java field), 16
width (Java field), 16
window (Java field), 16